


Liquid-Based Microbiology: the past, the present, the future

Matteo Miriani, Ph.D., Scientific Marketing Manager
 Copan Italia, Via Francesco Perotti, 10, 25125 Brescia BS

INTRODUCTION


In the last two decades, the world of microbiology underwent profound changes. Generally speaking, microbiology labs have seen an increased workload, and new pathogens have arisen; thus, new approaches needed to be developed and implemented, new diagnostics techniques have been developed, and a chronic lack of resources (both financial and human resources) forced labs to reshape and become more efficient and effective. All the phases of the diagnostics process were impacted in this process, including preanalytics.

Undoubtedly, one of the major improvements was the development of Liquid-Based Microbiology (LBM). This technology represented a departure from the conventional way to transport and preserve samples, from a solid format (a cotton or rayon swab inserted in a gelified transport media or a wet sponge) to a fully liquid solution, where the sample is released. Compared to the gel-based transport system, LBM introduces many benefits:



The uniform release of the sample into the transport media is one of the key features of LBM, and it is only possible by combining it with a FLOQSwabs®.


The benefits of moving from a gel-based transport system to LBM are not strictly related to the lab but can also be extended to the sample collection staff, the patients, and the procurement office, as noticeable in the table below:



All these improvements significantly impacted the microbiology diagnostics flow, making it easier, standardized, and automatable.

THE RIGHT PRODUCT FOR THE RIGHT APPLICATION

The modern microbiology lab is a mix of different technologies and platforms. Standardizing the sample collection using a few validated devices is possible and a real option. Said this, having the chance to choose different transport media characterized by different preservation specificities allows every professional to select the best transport media for any need.



This picture is meant only for flow demonstration purposes, and the data reported are not intended to replace the Instructions for Use (IFU). Always refer to the IFU for the final application compatibility or refer to a Copan representative. For molecular and antigen detection applications compatibility, refer to the device's manufacturer IFU.

CULTURE METHODS

The main objective for a bacteriology medium is to preserve the vitality of the microorganism that needs to be cultured, avoiding, at the same time, the overgrowth of commensal flora and high-fitness microorganisms. Intuitively, as the typology and physiology of the pathogens of human interest are very broad, designing a “universal media” is very complicated. Historically, at least three main gel-based transport systems were used according to the typology of bacteria to recover: Stuart or Amies gel w/o charcoal for aerobes, Stuart or Amies gel with charcoal for fastidious bacteria, and a specifically designed device for anaerobes. This high diversity of collection and transport devices definitively poses issues in terms of standardization and usability.

eSwab®

The Copan eSwab® was launched in 2006 and represented a fundamental improvement over the past transport system. In fact, with eSwab®, the microbiologist can standardize most of the needs on just one collection device. Copan eSwab® is

based on liquid Amies and is designed to preserve the vitality of aerobes bacteria, fastidious bacteria, and anaerobes for up to 48 hours at room temperature (apart from *N. gonorrhoea*, which is viable up to 24 hours). eSwab® was also proven to work efficiently for the viability preservation of fungi and yeast and for stabilizing DNA and antigens for molecular testing and antigen detection. eSwab® is the perfect combination with any automation platform for plate streaking and incubation.

eSwab® is currently used in routine by thousands of laboratories worldwide, and it was cited in over 180 peer-reviewed papers for multiple applications.

eSwab® was proven equivalent or superior to standard gel-based transport media in preserving bacteria viability and allowing Gram stain¹⁻². Results from stand-alone evaluations were excellent on aerobes on both mono species and poli species samples³⁻⁴, on fastidious bacteria⁵, and anaerobic bacteria⁶.

Antibiotic Susceptibility Testing (AST) is now among the most important – if not the most important – microbiology tests run in the lab. eSwab® fits this purpose perfectly, as it allows the collection of samples from different body sites, fully preserving the viability of the researched bacteria for potential antibiotic phenotypic detection. eSwab® was extensively evaluated, and it is currently routinely used worldwide for collecting samples to research MRSA from wounds⁷, skin, and respiratory sites⁸⁻⁹, as well as for the detection of CRE¹⁰. Also, eSwab® is considered the gold standard method by the CDC-Atlanta for collecting samples for *C. auris* and subsequent determination of antibiotic resistance patterns. Samples collected in eSwab® can be directly struck on selective plates and used as inoculum for enrichment broth.

eSwab® is the standard device for collecting *Streptococcus agalatie* (Strep B or GBS) in pregnant women¹¹⁻¹³ from samples from the anogenital area. The eSwab® ensures a long preservation time and easily allows the inoculum into selective pre-enrichment broth.

The eSwab® is also perfect for collecting and preserving samples for *Streptococcus pyogenes* (Strep A or GAS) to detect the leading cause of

bacterial pharyngitis and subsequent detection of antibiotic resistance¹⁴.

eSwab[®] has been proven extremely powerful in preserving the vitality of different types of fungi and yeast¹⁵. In this paper, the authors tested 19 isolates consisting of five yeast reference isolates (*Candida albicans* ATCC 10231, *Candida krusei* ATCC 6258, *Candida guilliermondii* ATCC 6260, *Candida glabrata* ATCC 66032, and *Cryptococcus neoformans* ATCC 66031), three reference dermatophytes (*Trichophyton mentagrophytes* ATCC 9533, *Trichophyton tonsurans* ATCC 28942, and *Trichophyton rubrum* ATCC 28188), five opportunistic hyaline molds (*Aspergillus niger*, *Aspergillus fumigatus*, *Lecytophthora sp.*, *Fusarium solani*, and *Trichosporon sp.*), three Zygomycetes (*Lichtheimia corymbifera*, *Mucor circinelloides*, and *Rhizopus microsporus*), and three threedematiaceous molds (*Curvularia clavata*, *Phialophora americana*, and *Alternaria alternata*). All the species were successfully recovered from the eSwab[®]. *Candida auris* was also successfully recovered and transported in eSwab[®] from the armpit and inguinal swabs.¹⁶

eSwab[®] is considered an excellent collection device for testing women's health syndromes and sexually transmitted diseases from vaginal sites in women, as well as from urethral sites from both women and male samples. *Trichomonas vaginalis* has been proven to be successfully recovered from eSwab[®]¹⁷. Among the different pathogens in this class, *Neisseria gonorrhoea* (NG) is one of the most concerning ones, especially because of the rise of antibiotic resistance in this pathogen. For a good epidemiological screening, the culturing of the pathogens is fundamental; thus, the vitality preservation in the collected sample is paramount. eSwab[®] was proven to preserve the vitality of NG^{17-19a} efficiently.

FecalSwab[®]

A critical diagnostic field for microbiology is the detection of gastrointestinal pathogens (GI). Infective gastroenteritis can be caused by different etiological causes: microorganisms, viruses, bacteria, and parasites can all cause acute or chronic diarrhea. In the case of bacterial pathogens detection in stool, if the analysis is

postponed after two hours from the sample collection (which is usually the case), the sample must be placed in a Cary-Blair-based preservative solution^{19b}. For this purpose, Copan developed the FecalSwab[®], a modified Cary-Blair medium specifically designed to preserve the vitality of entomopathogens at room temperature for up to 48 hours or up to 72 hours at 2-8 °C. The improvements of moving from native stool to FecalSwab[®] are multiple: higher sample stability, cleanliness, readiness of use, flexibility, and automatability. FecalSwab[®] can be used to preserve either native stool or rectal swabs and can be used for multiple serial or parallel analyses with different technologies.

Using FecalSwab[®] over an Amies or Stuart-based collection device is important because stool samples are extremely rich in commensal microorganisms, which may overgrow in a standard transport medium²⁰.

The use of FecalSwab[®] has been proven to double the *Campylobacter* recovery rate in clinical settings compared to native stool²¹. *Campylobacter* is probably one of the most fragile and faint bacteria causing enteritis, and the use of the FecalSwab[®] helps dramatically its recovery and, thus, diagnosis.

FecalSwab[®] was extensively evaluated also with automation platforms for plate streaking and incubation²², as well as with non-human stool veterinary samples²³⁻²⁴. Very interesting studies successfully evaluated FecalSwab[®] ability to isolate *Helicobacter spp* species from stool²⁵ and microbiome samples²⁶⁻²⁷.

Together with standard gastroenteritis diagnostics, the other most important field of application for the FecalSwab[®] is the preservation of rectal swabs for AMR screening, especially for ESBL, CRE, and VRE. The detection and monitoring of these resistances were evaluated by both cultural phenotypic assay and molecular approaches²⁸⁻³⁹.

C. difficile infection (CDI) is among the most important healthcare-acquired infections (HAI). The high mortality associated with CDI, the high risk of hospital outbreaks¹⁹, and the intrinsic complexity of *C. diff* diagnostics, stress, and even more, the need for an efficient preanalytical part.

FecalSwab® preserved stool was efficiently used to detect GDH and Toxin A and B⁴⁰⁻⁴¹.

UTM®

Nowadays, most of the diagnostics for viral infections are done by molecular assays or antigen tests. Despite this, the possibility of preserving the structure of the viruses remains paramount for allowing virus isolation for research, epidemiology, and clinical research. If the design of a transport media for preserving bacterial vitality recovery is complex, the creation of a proper viral transport media is even more critical as it must ensure that: 1) the viral capsid must stay intact; 2) the media must avoid bacterial and fungal growth 3) the media must be noncytotoxic to allow viral culture.

This is why only very few media can be considered “real” viral transport media, and Copan UTM can definitively be considered the gold standard. Over 200 publications used it for the recovery and later cultivation of multiple types of viruses, from respiratory to dermal and fecal.

MOLECULAR DIAGNOSTICS

The adoption of molecular biology by virology and microbiology labs is now a standard. In the early molecular days, these techniques required high skills and experience, dedicated areas for sample and reagent preparation, and experience in results interpretation; nowadays, with the latest automatic platform, molecular techniques are undoubtedly easier and more consistent than striking a plate. From small POC instruments to larger automated platforms that exploit different technologies (qPCR, LAMP, NASBA), nowadays, the offer is extremely high and differentiated. After the COVID-19 pandemic, even NGS and other sequencing techniques are becoming more present in diagnostics laboratories.

While for some pathogens, molecular assays are now practically the only choice (CT, NG, HPV...), for other pathogens, the option to go molecular or not can be driven by different factors like price, speed or reporting, or number of tests to run. Anyhow, the coexistence of molecular assays and classic microbiology techniques shall not be seen

as a dichotomy but as a combination able to improve modern laboratories’ results and supply complementary info. Also here, having the right sample quality is again paramount.

A good quality sample for molecular biology is when the nucleic acids (DNA and/or RNA) are preserved and not degraded. This may be achieved by preserving the cell's structure containing the nucleic acid, using active components that stabilize the DNA/RNA structure, or inhibiting the action of RNase and DNase.

In Copan’s LBM range, all the transport media are compatible with molecular biology assays, even with different specificity.

eSwab®

Copan eSwab® is validated to preserve the DNA for up to 5 days at room temperature. Even if the main intended use of eSwab® is preserving bacterial vitality, the flexibility of performing molecular biology on the same sample allows for standardized sample collection, with all the advantages arising from this.

eSwab® was successfully tested and validated for the detection of Strep A⁴², allowing the flexibility of molecular screening and culture reflex confirmation test and AST for the detection of Group B streptococcus⁴³, which allows standardizing on just one collection device both flow of patients that undergo the screening between the 35 to the 37 weeks of pregnancy by culture, and the ones tested in emergency during the labor by molecular.

eSwab® was successfully used with different molecular assays for the detection of *Chlamydia trachomatis* and of *Neisseria gonorrhoea* (NG)⁴⁵⁻⁴⁶, allowing also, in this case, the subsequential NG culture for AST. eSwab® can be used for the collection of samples for the molecular detection of women's health syndromes, like vaginosis and vaginitis syndromic panels⁴⁴.

As mentioned in the previous chapter on using eSwab® for resistant pathogens by culture, the same approach can be done by molecular assays. eSwab® is currently used to detect the genetic determinants of MRSA, CRE, and VRE⁴⁵⁻⁴⁷. Also, in this case, eSwab® allows the standardization of the collection part while leaving the freedom to

adapt the detection technology according to the case.

FecalSwab®

One of the fastest-growing fields where molecular biology plays a more critical role is the detection of gastrointestinal pathogens. MDx allows fast and comprehensive screening of multiple pathogens, including viruses, bacteria, and parasites. This is very important to preserve the vitality of the bacteria eventually present in the sample and to allow the culture and isolation of positive samples for bacteria.

In this landscape, FecalSwab® plays a crucial role. FecalSwab® is perfect for collecting stool samples for PCR screening of different bacteria, viruses, and parasites with different multiplex assays⁴⁸⁻⁵⁵. Typically, *C. difficile* diagnostics is managed out of the standard GI flow, and FecalSwab® has been proven to be highly efficient for collecting samples and preserving stools and compatible with the primary molecular assays available on the market⁵⁶⁻⁶⁰.

FecalSwab® is extremely well-fitted to preserve samples intended for genetic determinants of antibiotic resistances, especially ESBL, CRE, and VRE from rectal swab²⁸⁻³⁹.

A final mention of the FecalSwab® molecular capabilities is the detection of gastroenteritis pathogens from rectal swabs instead of stool⁶¹⁻⁶³. Different studies indeed demonstrated that the performances of the two types of samples are comparable. Using rectal swabs is highly advantageous for specific populations (newborns, children, and elderly) where stool collection can be problematic. In contrast, the collection of a rectal swab is fast and immediate, allowing prompt results delivery for a faster treatment response.

UTM®

When clinical necessity is the detection of viruses, the most known product is UTM®. Even if the UTM®'s original intended use is to preserve the viral integrity, most of the samples collected in UTM® are analyzed by NAAT or by antigen assay. UTM® is routinely used for collecting samples for the detection of respiratory viruses (from either

the upper or the lower respiratory tract), STI viruses (HSV 1&2 from blister and sores, monkeypox virus), and bacteria (*Chlamydia spp.*, *Mycoplasma spp.*, *Ureoplasma spp.*), conjunctivitis viruses, and stool samples.

CONCLUSIONS

Liquid-Based Microbiology, together with the FLOQSwabs®, represented the most important improvement in the microbiology preanalytical phase of the last 100 years. The Copan's LBM line has now been used by thousands of laboratories worldwide, in millions of pieces per year. The quality and flexibility of the Copan's LBM allow one to find the right product for the application or pathogen you need to detect.

NEED MORE INFO?

Visit our web-site <https://www.copangroup.com/> or contact us at info@copangroup.com.

REFERENCES

- 1) Comparison of Copan eSwab with the Copan Venturi Transystem for the quantitative survival of *Escherichia coli*, *Streptococcus agalactiae* and *Candida albicans* S Nys, S Vijgen, K Magerman, R Cartuyvels Eur J Clin Microbiol Infect Dis . 2010 Apr;29(4):453-6.
- 2) Comparison of the eSwab collection and transportation system to an amies gel transystem for Gram stain of clinical specimens Carla Fontana, Marco Favaro, Dolores Limongi, Jana Pivonkova, Cartesio Favalli BMC Res Notes. 2009; 2: 244.
- 3) Comparison of the Copan ESwab system with two Amies agar swab transport systems for maintenance of microorganism viability Kenneth G Van Horn, Carol D Audette, Denise Sebeck, Kelly A Tucker J Clin Microbiol. 2008 May;46(5):1655-8. doi: 10.1128/JCM.02047-07.
- 4) Evaluation of three swab transport systems for the maintenance of clinically important bacteria in simulated mono- and polymicrobial samples Eva Tano, Asa Melhus APMIS . 2011 Mar;119(3):198-203. doi: 10.1111/j.1600-0463.2010.02710.x.
- 5) Comparison of the Copan eSwab System with an Agar Swab Transport System for Maintenance of Fastidious Anaerobic Bacterium Viability Kerin L Tyrrell , Diane M Citron, Eliza S Leoncio, Ellie J C Goldstein J Clin Microbiol . 2016 May;54(5):1364-7. doi: 10.1128/JCM.03246-15.
- 6) Extensive evaluation of fastidious anaerobic bacteria recovery from the Copan eSwab® transport system Thomas Demuyser, Deborah De Geyter, Daisy Van Dorpe, Kristof Vandoorslaer, Ingrid Wybo J Microbiol Methods . 2018 Jan;144:73-78. doi: 10.1016/j.mimet.2017.11.009. Epub 2017 Nov 13.
- 7) Clinical evaluation of the Copan eSwab for methicillin-resistant *Staphylococcus aureus* detection and culture of wounds V Saegeman, J Flamaing, J Muller, W E Peetermans, J Stuyck, J Verhaegen Eur J Clin Microbiol Infect Dis . 2011 Aug;30(8):943-9. doi: 10.1007/s10096-011-1178-1.
- 8) Evaluation of the Copan eSwab transport system for the detection of methicillin-resistant *Staphylococcus aureus*: a laboratory and clinical study Annick Smismans, Jan Verhaegen, Annette Schuermans, Johan Frans Diagn Microbiol Infect Dis . 2009 Oct;65(2):108-11.
- 9) Nasal Screening for MRSA: Different Swabs – Different Results! Philipp Warnke, Hagen Frickmann, Peter Ottl, Andreas Podbielski PLoS One. 2014; 9(10): e111627.
- 10) Laboratory evaluation of the ESwab transport system for the recovery of carbapenem-resistant *Acinetobacter baumannii* J Moran-Gilad, D Schwartz, S Navon-Venezia, Y Carmeli Eur J Clin Microbiol Infect Dis . 2012 Jul;31(7):1429-33. doi: 10.1007/s10096-011-1460-2.
- 11) Evaluation of optimal storage temperature, time, and transport medium for detection of group B *Streptococcus* in StrepB carrot broth Ashton Trotman-Grant, Trisha Raney, Jennifer Dien Bard J Clin Microbiol . 2012 Jul;50(7):2446-9. doi: 10.1128/JCM.00238-12.
- 12) Clinical evaluation of the walk-away specimen processor and ESwab for recovery of *Streptococcus agalactiae* isolates in prenatal screening specimens Blake W Buchan, Wendy J Olson, Tami-Lea A Mackey, Nathan A Ledebor Clin Microbiol. 2014 Jun;52(6):2166-8. doi: 10.1128/JCM.00374-14.
- 13) Comparison of different sampling techniques and of different culture methods for detection of group B streptococcus carriage in pregnant women Nabil A El Aila, Inge Tency, Geert Claeys, Bart Saerens, Piet Cools, Hans Verstraelen, Marleen Temmerman, Rita Verhelst, Mario Vaneechoutte BMC Infect Dis . 2010 Sep 29;10:285. doi: 10.1186/1471-2334-10-285
- 14) The Utility of Rapid Group A *Streptococcus* Molecular Testing Compared with Throat Culture for the Diagnosis of Group A Streptococcal Pharyngitis in a High-Incidence Rheumatic Fever Population Amanda Taylor, Susan Morpeth, Rachel Webb, Susan Taylor J Clin Microbiol . 2021 Nov 18;59(12):e0097821. doi: 10.1128/JCM.00978-21.
- 15) Evaluation of the Copan ESwab Transport System for Viability of Pathogenic Fungi by Use of a Modification of Clinical and Laboratory Standards Institute Document M40-A2 Bharat Gandhi, Richard Summerbell, Tony Mazzulli J Clin Microbiol . 2018 Jan 24;56(2):e01481-17. doi: 10.1128/JCM.01481-17
- 16) Identification of Drug Resistant *Candida auris* Milena Kordalewska and David S. Perlin Front Microbiol. 2019; 10: 1918. Published online 2019 Aug 20. doi: 10.3389/fmicb.2019.01918
- 17) Viability of *Trichomonas vaginalis* in Copan Universal Transport Medium and eSwab Transport Medium C. A. Rivers, J. R. Schwebke J Clin Microbiol. 2008 Sep; 46(9): 3134–3135.
- 18) Evaluation of the Copan eSwab®, a liquid-based microbiology transport system, for the preservation of *Neisseria gonorrhoeae* at different temperatures Lindy Gumede , Frans Radebe , Duduzile Nhlapo , Venessa Maseko , Tendesayi Kufa-Chakezha and Ranmini Kularatne Southern African Journal of Infectious Diseases 2017; 32(3):96–99 <https://doi.org/10.1080/23120053.2017.1313935>
- 19a) Antimicrobial susceptibility of *Neisseria gonorrhoeae* isolates and syndromic treatment of men with urethral discharge in Kingston, Jamaica, 2018–19 Suzette M Cameron-McDermott, Geoffrey J Barrow, Alicia M Webster, Carrington O De La Haye, Denise H E Wood, Violet M Lewis, Alison Nicholson, Glendee Y Reynolds-Campbell, Camille-Ann A Thoms-Rodriguez, Karen J Roye-Green, Nakeisha Otto-Stewart, Zahra N Miller, Jennifer A Tomlinson , Nicola Skyers, Magnus Unemo, Joshua J Anzinger J Antimicrob Chemother . 2021 Dec 24;77(1):218-222. doi: 10.1093/jac/dkab340.
- 19b) Laboratory diagnosis of bacterial gastroenteritis Romney M Humphries 1, Andrea J Linscott Clin Microbiol Rev. 2015 Jan;28(1):3-31. doi: 10.1128/CMR.00073-14.
- 20) Comparison of fecalswab and eswab devices for storage and transportation of diarrheagenic bacteria Jari J Hirvonen, Suvi-Sirkku Kaukoranta J Clin Microbiol . 2014 Jul;52(7):2334-9. doi: 10.1128/JCM.00539-14.
- 21) Enhanced culture recovery of campylobacter with modified cary-blair medium: a practical field experience Clémence Massip, Hélène Guet-Revillet, Marion Grare, Agnès Sommet, Damien Dubois J Microbiol Methods. 2018 Jun;149:53-54. doi: 10.1016/j.mimet.2018.05.001.
- 22) Evaluating the preservation and isolation of stool pathogens using the copan fecalswab™ transport system and walk-away specimen processor Lee W Goneau, Anthony Mazzulli, Xhermi Trimi, Ana Cabrera, Pauline Lo, Tony Mazzulli Diagn Microbiol Infect Dis . 2019 May;94(1):15-21. doi: 10.1016/j.diagmicrobio.2018.11.020
- 23) Occurrence of salmonella enterica in grey-headed flying foxes from New South Wales F McDougall, M Power Australian Veterinary Journal, 2021 <https://doi.org/10.1111/avj.13116>
- 24) Detection, Identification, and Antimicrobial Susceptibility of *Campylobacter* spp. and *Salmonella* spp. from Free-ranging Nonhuman Primates in Sri Lanka Cecilia Tegner, N. P. Sunil-Chandra, W. R. P. L. I. Wijesooriya, B. Vijitha Perera, Ingrid Hansson, Åsa Fahlman J Wildl Dis (2019) 55 (4): 879–884. <https://doi.org/10.7589/2018-08-199>
- 25) *Helicobacter cinaedi* bacteraemia secondary to enterocolitis in an immunocompetent patient Sofie Larsen Rasmussen, Iben Ørsted, Irene Harder Tarpgaard, Hans Linde Nielsen Gut Pathog . 2021 Apr 23;13(1):26. doi: 10.1186/s13099-021-00422-8.
- 26) The effect of sampling and storage on the fecal microbiota composition in healthy and diseased subjects Danyta I Tedjo, Daisy M A E Jonkers, Paul H Savelkoul, Ad A Masclee, Niels van Best, Marieke J Pierik, John Penders PLoS One. 2015 May 29;10(5):e0126685. doi: 10.1371/journal.pone.0126685.
- 27) Usability of rectal swabs for microbiome sampling in a cohort study of hematological and oncological patients Lena M. Biehl, Debora Garzetti, Fedja Farowski, Diana Ring, Martin B. Koepfel, Holger Rohde, Philippe Schafhausen, Bärbel Stecher, Maria J. G. T.

Vehreschild PLoS One April 15, 2019
<https://doi.org/10.1371/journal.pone.0215428>

28) Rapid and accurate exdr screening: use xpert carba-r® with FecalSwab® Eric Farfour, Alexandra Lomont, Vincent Fihman, Marion Lecuru, Sophie Hüssler, Souad Ouzani, Jean-Ralph Zahar, Laurent Dortet *Diagn Microbiol Infect Dis* . 2021 Apr;99(4):115279. doi: 10.1016/j.diagmicrobio.2020.115279.

29) Multi-drug-resistant enterobacteriaceae carriage in highly exposed nursing homes: prevalence in western France G. Birganda, N. Hayatgheibc, P. Bemerc, V. Guilloteau, C. Legeay. S. Perron, G. Chapelet, S. Corvece, C. Bourigault, E. Batardec, D. Lepelletier *Journal of hospital infection*, volume 97, issue 3, p258-259, november 2017, volume 97, issue 3, p258-259, November 2017

30) High prevalence of multidrug-resistant Enterobacteriaceae among residents of long term care facilities in Amsterdam, the Netherlands Eline van Dulm, Aletta T. R. Tholen, Annika Pettersson, Martijn S. van Rooijen, Ina Willemsen, Peter Molenaar, Marjolein Damen, Paul Gruteke, Paul Oostvogel, Ed J. Kuijper, Cees M. P. M. Hertogh, Christina M. J. E. Vandenbroucke-Grauls, Maarten Scholing *PLoS One*, 2019, 14.9: e0222200.
<https://doi.org/10.1371/journal.pone.0222200>

31) Genomic and phenotypic analysis of an esbl-producing e. coli st1159 clonal lineage from wild birds in Mongolia Peter Schierack, Stefan E. Heiden, , Muhammad Moman Khan, Lena Nikolaus, Rafal Kolenda, Michael Stubbe, Davaa Lkhagvasuren, Stefan Rödiger, Sebastian Guenther, Katharina Schaufler *Front. Microbiol.*, 21 July 2020, Sec. Antimicrobials, Resistance and Chemotherapy, Volume 11 - 2020 | <https://doi.org/10.3389/fmicb.2020.01699>

32) Gastrointestinal carriage of antimicrobial resistance in school-aged children in three municipalities of timor-leste Tessa Oakley, Brandon Le, Virginia da Conceicao, Ian Marr, Carolina Maia , Messias Soares, Joana Correia Belo, Nevio Sarmento, Endang da Silva , Salvador Amaral, Susana Vaz Nery, Sarah Lynar, Joshua R. Francis, Jennifer Yan *Antibiotics* 2022, 11(9), 1262; <https://doi.org/10.3390/antibiotics11091262>

33) Risk factors for community colonization with extended-spectrum cephalosporin-resistant enterobacterales (escr) in botswana: an antibiotic resistance in communities and hospitals (arch) study Ebbing Lautenbach, Mosepele Mosepele, Rachel M Smith, Ashley Styczynski, Robert Gross, Leigh Cressman, Anne Jaskowiak-Barr, Kevin Alby, Laurel Glaser, Melissa Richard-Greenblatt, Laura Cowden, Kgotlaetsile Sewawa, Dimpho Otukile, Giacomo M Paganotti, Margaret Mokomane, Warren B Bilker, Naledi Mannathoko *Clinical Infectious Diseases*, Volume 77, Issue Supplement_1, 1 July 2023, Pages S89–S96, <https://doi.org/10.1093/cid/ciad259>

34) Performance of the bd max™ instrument with check-direct cpe real-time pcr for the detection of carbapenemase genes from rectal swabs, in a setting with endemic dissemination of carbapenemase-producing Enterobacteriaceae Alberto Antonelli, Fabio Arena, Tommaso Giani, Olga Lorenza Colavecchio, Stoyanka Valentinova Valeva, Suzanne Paule, Peter Boleij, Gian Maria Rossolini *Diagnostic Microbiology and Infectious Disease* Volume 86, Issue 1, September 2016, Pages 30-34 <https://doi.org/10.1016/j.diagmicrobio.2016.06.002>

35) Longitudinal monitoring reveals the emergence and spread of blaGES-5-harboring carbapenem-resistant Klebsiella quasipneumoniae in a Hong Kong hospital wastewater discharge line Xin Liu a b c, Matthew K.L. Wong, Dengwei Zhang, Darren C.L. Chan, Olivia S.K. Chan, Gary P.L. Chan, Marcus Ho-Hin Shum, Ye Peng, Christopher K.C. Lai, Benjamin J. Cowling, Tong Zhang, Keiji Fukuda, Tommy Tsam-Yuk Lam, Hein M. Tun *Science of The Total Environment* Volume 903, 10 December 2023, 166255 <https://doi.org/10.1016/j.scitotenv.2023.166255>

36) Colonization with extended-spectrum cephalosporin-resistant enterobacterales (escr) and carbapenem-resistant enterobacterales (cre) in healthcare and community settings in botswana: an antibiotic

resistance in communities and hospitals (arch) study Naledi Mannathoko, Mosepele Mosepele, Robert Gross, Rachel M. Smith, Kevin Alby, Laurel Glaser, Melissa Richard-Greenblatt, Rebekah Dumm, Aditya Sharma, Anne Jaskowiak-Barr, Leigh Cressman, Kgotlaetsile Sewawa, Laura Cowden, Emily Reese, Dimpho Otukile, Giacomo M. Paganotti, Margaret Mokomane, Ebbing Lautenbach *International Journal of Infectious Diseases* Volume 122, September 2022, Pages 313-320 <https://doi.org/10.1016/j.ijid.2022.06.004>

37) Destination shapes antibiotic resistance gene acquisitions, abundance increases, and diversity changes in Dutch travelers Alaric W. D'Souza, Manish Boolchandani, Sanket Patel, Gianluca Galazzo, Jarne M. van Hattem, Maris S. Arcilla, Damian C. Melles, Menno D. de Jong, Constance Schultsz, COMBAT Consortium, Gautam Dantas , John Penders *Genome Med* 13, 79 (2021). <https://doi.org/10.1186/s13073-021-00893-z>

38) Outbreak of ceftazidime-avibactam resistant Klebsiella pneumoniae carbapenemase (KPC)-producing Klebsiella pneumoniae in a COVID-19 intensive care unit, Italy: urgent need for updated diagnostic protocols of surveillance cultures G. Bianco, M. Boattini, A. Bondi, S. Comini, T. Zaccaria, R. Cavallo, C. Costa *Journal of Hospital Infection*, VOLUME 122, P217-219, APRIL 2022, DOI:<https://doi.org/10.1016/j.jhin.2022.02.001>

39) Implementation of chromatic super caz/avi® medium for active surveillance of ceftazidime-avibactam resistance: preventing the loop from becoming a spiral Gabriele Bianco, Matteo Boattini, Sara Comini, Alessio Leone, Alessandro Bondi, Teresa Zaccaria, Rossana Cavallo, & Cristina Costa *European Journal of Clinical Microbiology & Infectious Diseases*, 41(9), 1165-1171.

40) Evaluation of liaison® c. difficile glutamate dehydrogenase and liaison® c. difficile toxin a and b in copan fecalswab™ samples in a three-step algorithm for the diagnosis of Clostridium difficile infection Massimo Oggioni, Alessandra Bielli, Alice Nava and Daniela Adele Pia Campisi *Microbiologia Medica* June 201530(1) DOI:10.4081/mm.2015.4953

41) Prospective Evaluation of the mariPOC Test for Detection of Clostridioides difficile Glutamate Dehydrogenase and Toxins A/B Roosa Savolainen, Juha M. Koskinen, Silja Mentula, Janne O. Koskinen, Suvi-Sirkku Kaukoranta *Journal of Clinical Microbiology* Vol. 58, No. 4 <https://doi.org/10.1128/jcm.01872-19>

42) Comparison of illumigene Group A Streptococcus Assay with Culture of Throat Swabs from Children with Sore Throats in the New Zealand School-Based Rheumatic Fever Prevention Program Arlo Upton, Liselle Bissessor, Elizabeth Farrell, Stanford T Shulman, Xiaotian Zheng, Diana Lennon *J Clin Microbiol* 2016 Jan;54(1):153-6. doi: 10.1128/JCM.02440-15.

43) Detection of Group B Streptococcus Directly from Collected ESwab Samples by Use of the BD Max GBS Assay Suzane Silbert, Talita T Rocchetti, Alicia Gostnell , Carly Kubasek, Raymond Widen *J Clin Microbiol* . 2016 Jun;54(6):1660-1663. doi: 10.1128/JCM.00445-16.

44) Diagnostic performance of two molecular assays for the detection of vaginitis in symptomatic women Alexandra Thompson, Karen Timm, Noelle Borders, Liz Montoya, Karissa Culbreath *Eur J Clin Microbiol Infect Dis* . 2020 Jan;39(1):39-44.

45) Successful Combination of Nucleic Acid Amplification Test Diagnostics and Targeted Deferred Neisseria gonorrhoeae Culture Carolien M. Wind, Henry J. C. de Vries, Maarten F. Schim van der Loeff, Magnus Unemo, e and Alje P. van Dam *J Clin Microbiol*. 2015 Jun; 53(6): 1884–1890.

46) Copan ESwab™, a Liquid Based Microbiology Device, Can Be Used For the Preservation of Neisseria Gonorrhoeae For Culture and For Detection of CT/NG by Gen-Probe® APTIMA® COMBO 2 Assay M Jett-Goheen, M Barnes, D Duncan, N Quinn, K Ghanem, C Gaydos <http://dx.doi.org/10.1136/sextrans-2013-051184.1113>

45) Rapid identification from rectal swabs of the clinically most relevant carbapenemase genes from gram-negative bacteria using the BD MAX Check-Points CPO Assay Sergio García-Fernández, Patricia J Simner, Gina Thomson, Matthew Faron, Ron Bosboom,

- Arjanne van Griethuijsen, María García-Castillo, Renee Harris, Nathan A Ledebøer, Rafael Cantón, Kenneth S Thomson *Diagn Microbiol Infect Dis.* 2022 Jan;102(1):115554.
- 46) Detection of colonization by carbapenem-resistant organisms by real-time polymerase chain reaction from rectal swabs in patients with chronic renal disease T F T Rezende, A M Doi, M G Quiles, A C C Pignatari, S Manfrendi, C Grothe, M Taminato, D A Barbosa *J Hosp Infect.* 2017 Jun;96(2):123-128. doi: 10.1016/j.jhin.2017.03.021.
- 47) Development of a Real-Time PCR Protocol Requiring Minimal Handling for Detection of Vancomycin-Resistant Enterococci with the Fully Automated BD Max System Alexander H Dalpke, Marjeta Hofko, Stefan Zimmermann *J Clin Microbiol.* 2016 Sep;54(9):2321-9. doi: 10.1128/JCM.00768-16.
- 48) Evaluation of copan fecalswab™ preserved stool specimens with the bd max™ enteric bacterial panel and the bd max™ extended enteric bacterial panel Hanen Fernandez Rojas, Amorce Lima, Carly Kubasek, Alicia Gostnell, Suzane Silbert *Diagn Microbiol Infect Dis.* 2020 Aug;97(4):115055
- 49) Evaluation of the new fecalswab system for maintaining stability of stool samples submitted for molecular tests Suzane Silbert, Alicia Gostnell, Carly Kubasek, Raymond Widen *J Clin Microbiol.* 2017 May;55(5):1588-1590. doi: 10.1128/JCM.00273-17.
- 50) Evaluation of the fecalswab for stool specimen storage and molecular detection of enteropathogens on the bd max system Melissa Richard-Greenblatt, Candy Rutherford, Kathy Luinstra, Ana María Cárdenas, Xiaoli Lilly Pang, Padman Jayaratne, Marek Smieja *J Clin Microbiol.* 2020 Aug 24;58(9):e00178-20. doi: 10.1128/JCM.00178-20
- 51) Comparison of a Rapid Multiplex Gastrointestinal Panel with Standard Laboratory Testing in the Management of Children with Hematochezia in a Pediatric Emergency Department: Randomized Controlled Trial Jianling Xie, Kelly Kim, Byron M. Berenger, Linda, Otto G. Vanderkooi, Silviu Grisaru, Stephen B. Freedman *Microbiology Spectrum* (2023): e00268-23. <https://doi.org/10.1128/spectrum.00268-23>
- 52) Evaluation of the Allplex™ GI-Helminth (I) Assay, the first marketed multiplex PCR for helminth diagnosis. Autier, B., Gangneux, J. P., & Robert-Gangneux, F. *Parasite*, 28. (2021). doi: 10.1051/parasite/2021034
- 53) The clinical significance of *dientamoeba fragilis* and *blastocystis* in human stool—retrospective cohort study Shasha, D., Grupel, D., Treigerman, O., Prajgrod, G., Paran, Y., Hacham, D., Zacay, G. *Clinical Microbiology and Infection.* 2023 <https://doi.org/10.1016/j.cmi.2023.09.003>
- 54) Carriage of *Blastocystis* spp. in travellers - A prospective longitudinal study van Hattem, J. M., Arcilla, M. S., Schultsz, C., Bootsma, M. C., Verhaar, N., Rebers, S. P., ... & Verbrugh, H. A. *Travel Medicine and Infectious Disease*, 27, 87-91. <https://doi.org/10.1016/j.tmaid.2018.06.005>
- 55) Easyscreen™ Enteric Protozoa Assay For the Detection of Intestinal Parasites: A Retrospective Bi-Center Study Giorgio Dirani; Silvia Zannoli; Elena Paesini; Patrizia Farabegoli; Barbara Dalmo; Caterina Vocale; Giovanna Liguori; Stefania Varani; Vittorio Sambri *J Parasitol* (2019) 105 (1): 58–63. <https://doi.org/10.1645/18-52>
- 56) Direct detection and characterization of *clostridium difficile* from a novel collection device to improve laboratory workflow Fabrizio Barbanti, Patrizia Spigaglia *APMIS.* 2019 Jun;127(6):449-453. doi: 10.1111/apm.12942.
- 57) Performances of the bd max™ cdiff assay for the detection of toxigenic *clostridioides difficile* using cary-blair preserved samples Ranc Anne-Gaëlle, Sentis Célie, Couturier Jeanne, Barbut Frédéric, Tristan Anne, Buis Coralie, Santos Kevin, Farrah Sonia, Duraffour Pascaline, Vandenesch Francois, Dauwalder Olivier, Laurent Frédéric *Diagn Microbiol Infect Dis.* 2022 Jul;103(3):115701. doi:10.1016/j.diagmicrobio.2022.115701.
- 58) Novel portable platform for molecular detection of toxigenic *Clostridium difficile* in faeces: a diagnostic accuracy study Hirvonen, J. J., Matero, P., Siebert, C., Kauppila, J., Vuento, R., Tuokko, H., & Boisset, S. *European Journal of Clinical Microbiology & Infectious Diseases*, 36, 783-789. (2017).
- 59) Comparison of bd max cdiff and genomera c. *difficile* molecular assays for detection of toxigenic *clostridium difficile* from stools in conventional sample containers and in fecalswabs J J Hirvonen, S-S Kaukoranta *Eur J Clin Microbiol Infect Dis.* 2015 May;34(5):1005-9. doi: 10.1007/s10096-015-2320-2.
- 60) Evaluation of copan fecalswab as specimen type for use in xpert c. *difficile* assay Michael J Mashock, Matthew L Faron, Blake W Buchan, Nathan A Ledebøer *J Clin Microbiol.* 2017 Oct;55(10):3123-3129. doi: 10.1128/JCM.00369-17.
- 61) enteropathogen detection in children with diarrhoea, or vomiting, or both, comparing rectal flocced swabs with stool specimens: an outpatient cohort study Stephen B Freedman, Jianling Xie, Alberto Nettel-Aguirre, Bonita Lee, Linda Chui, Xiao-Li Pang, Ran Zhuo, Brendon Parsons, James A Dickinson, Otto G Vanderkooi, Samina Ali, Lara Osterreicher, Karen Lowerison, Phillip I Tarr *Lancet Gastroenterol Hepatol.* 2017 Sep;2(9):662-669. doi: 10.1016/S2468-1253(17)30160-7
- 62) Evaluation of rectal swab use for the determination of enteric pathogens: a prospective study of diarrhea in adults T Kotar, M Pirš, A Steyer, T Cerar, B Šoba, M Skvarc, M Poljšak Prijatelj, T Lejko Zupanc *Clin Microbiol Infect.* 2019 Jun;25(6):733-738. doi: 10.1016/j.cmi.2018.09.026.
- 63) Rotavirus detection in bulk stool and rectal swab specimens in children with acute gastroenteritis in Norway Moustafa Gibory, Ildri Haltbakk, Elmira Flem, Kirsti Vainio, Beatriz Valcarcel Salamanca, Ketil Størdal, Svein Arne Nordbø, Kirsti Jakobsen, Elisebet Haarr, Susanne Gjeruldsen *Dudman Journal of Clinical Virology* Volume 97, December 2017, Pages 50-53 <https://doi.org/10.1016/j.jcv.2017.10.017>