



# Whitley Modified Atmosphere Systems

## Product Comparison Chart

			
	<b>Instant Access Porthole System</b>	<b>Manual Porthole System 1</b>	<b>Manual Porthole System 2</b>
<b>Workstation Model</b>	<ul style="list-style-type: none"> <li>Eliminates sleeves and cuffs</li> <li>Instant access</li> <li>International patents</li> <li>Ideal for more often / shorter duration use</li> </ul>	<ul style="list-style-type: none"> <li>20 Sample transfer capacity</li> <li>Sleeve and cuff attachment</li> <li>Oval, ergonomic shape</li> <li>Integrated shelf for sample transfer</li> </ul>	<ul style="list-style-type: none"> <li>Sleeve and cuff attachment</li> <li>Oval, ergonomic shape</li> <li>Ideal for less often / longer duration use</li> </ul>
A20	—	●	—
A25	○	—	○
A35	○	—	○
A45	▲ ○	—	▲ ○
A55	○	—	○
A85	—	●	—
A95	—	●	—
A135	—	—	●
A155	—	—	●
M35	○	—	○
M45	—	—	●
M85	—	●	—
M95	—	●	—
H35	—	—	●
H45	—	—	●
H85	—	●	—
H95	—	●	—
H135	—	—	●
H155	—	—	●
i2	—	—	●

**KEY:**

— Not available   ● Fitted as standard   ○ Can be specified   ▲ Far left hand port always Instant Access