

IMPORTANT INFORMATION TO READ and RETURN

Installation Requirements for Protos 3

Thank you for choosing one of our products for your laboratory. To enable our engineers to perform an efficient, trouble-free installation please study, complete and fax this form to us on **01274 531197**. Should you have any questions, please do not hesitate to contact us, as we are here to help. When we have received the completed form, our Service Department will contact you to arrange a mutually convenient installation date.

**The following information represents the ideal requirement.
Please contact us IMMEDIATELY if your intended location does not match this specification.**



Space Requirements

The weight of the equipment is 20kg. The bench allocated must be flat, level and of sufficient size to support the base fully.

Bench Space Required

Depth	Length	Height
<i>mm</i>	<i>mm</i>	<i>mm</i>
500	400	540

Minimum PC Requirements

- Processor: Intel Core 2 or Core i processor 2.13 GHz or better
- Memory: 2GB (4GB recommended)
- Hard drive: 160GB (or larger)
- Graphics: Intel X4500 or better (Windows 7 Experience rating for graphics should be a minimum of 3.0, ideally higher)
- Operating system: Windows XP Professional 32 bit version Service Pack 3 or Windows 7 Professional 32 bit version or Windows 8.1 Professional. (Home versions not suitable. Windows Vista not supported).
- Two available high-speed USB 2.0 ports (one for camera, one for peripherals hub)
- Monitor: 1280 x 1024 Minimum resolution (not compatible with lower resolutions as some elements of the application become inaccessible)

Note: The combination of the minimum display resolution and the USB ports will preclude most laptops.

Mains Requirements

Electricity Supply
230 +/- 10% V

Sockets
2 x Three Pin, 13 amp

Delivery and installation are free of charge (unless otherwise agreed). If our engineers are unable to install the unit and a return journey is necessary, **a charge may be made.**

It is essential that this form is completed and returned, to avoid delay to your installation.

THANK YOU FOR THINKING WHITLEY

Signature

Title

Print Name

Establishment