

Whitley Workstations

Single Plate Entry Systems

Features:

Introduce individual samples without compromising internal atmospheric conditions

Versions for the majority of our workstations can be retro fitted

Perfect for on-call and emergency use



The Whitley Workstation range includes a variety of options and accessories enabling you to tailor a system to your needs and your laboratory working practices.

A very useful option is a single plate entry system, which allows individual Petri dishes to be introduced into the chamber without compromising internal atmospheric conditions. Multiwell plates, cell culture flasks and similar sized objects can be introduced using the same system.

The single plate entry system can be extremely useful for on-call and emergency use where individual samples can be introduced quickly without the need to use the airlock or ports.

More convenient working methods and lower gas usage, what more could you want?

	9cm				15cm
	A03125	A03225	A03325	A07208	A06103
Workstation Model	w: 107mm, h: 26mm Front loading Can be retro fitted	w: 107mm, h: 26mm Front loading Longer tray for larger items Can be retro fitted	w: 98mm, h: 20mm Front loading With a folding tray for an unobstructed work area Can be retro fitted	w: 107mm, h: 26mm Front loading Can be retro fitted on units with a removable front	w: 152mm, h: 52mm Located on the left hand panel of the unit Can be retro fitted
A20	-	-	-	-	●
A25	-	-	-	-	●
A35	-	-	-	●	●
A45	-	-	-	●	●
A55	-	-	-	●	-
A85	●	●	●	-	-
A95	●	●	●	-	-
A135 H	-	-	-	●	●
A135 GMP	-	-	-	●	●
A155 H	-	-	-	●	●
A155 GMP	-	-	-	-	-
M35	-	-	-	●	●
M45	-	-	-	●	●
M85	●	●	●	-	-
M95	●	●	●	-	-
M135 H	-	-	-	●	●
M155 H	-	-	-	●	●
H25	-	-	-	-	●
H35	-	-	-	●	●
H45	-	-	-	●	●
H85	●	●	●	-	●
H95	●	●	●	-	●
H135 H	-	-	-	●	●
H135 GMP	-	-	-	-	-
H155 H	-	-	-	●	-
H155 GMP	-	-	-	-	-
i2	-	-	-	-	-
KEY:	- Not available			● Option available	

Single Plate Entry Systems are designed for occasional use and not as a direct alternative to airlocks or porthole systems that allow sample transfer. *accommodates 90mm Petri dishes & items of a similar size.

In the interests of a policy of continuous product improvement the company reserves the right to alter product specifications without prior notice. All rights reserved. © 2023 Don Whitley Scientific Limited.